

STEP 1

INGREDIENTS		PREPARATION
JOYBASE GRANITA NEUTRA	g 1.000	Procedure fro Granita with slush machine
WATER	g 6.000-7.000	Mix all the ingredients until obtaining a smooth consistency. Pour in the slush machine and start
JOYPASTE IN THE CHOSEN FRUIT	g 100/kg acqua	the cooling system.
		Procedure for Granita in tub

Mix all the ingredients until obtaining a smooth consistency. Pour in the batch freezer and start the appropriate program for Granita (if not available, remove the product while it is still liquid).

Pour the Granita in the tub and stir often.

